




INTRODUCTION

PCT (procalcitonin) is a potential biomarker of sepsis. Serum levels increase during bacterial infections allowing differentiation from non-infectious diseases.

The dog PCT assay uses SPARCL™¹ technology. Two PCT antibody conjugates are used; one to horseradish peroxidase (HRP), the other to acridan, a chemiluminescent substrate. When HRP and acridan conjugates bind to PCT they are brought into proximity. Upon addition of hydrogen peroxide, HRP catalyzes oxidation of proximal acridan molecules causing a flash of luminescence that is proportional to PCT concentration.

Diluted serum samples and standards (100 µl) are dispensed into test tubes and mixed with 0.5 ml of combined acridan and HRP conjugates. After 45 minutes, tubes are placed in the VetBio-1 luminometer. Luminescence is measured after injection of background reducer (50 µl), that eliminates nonspecific luminescence, followed by trigger-solution containing hydrogen peroxide (100 µl). The concentration of PCT is determined from the ratio of sample luminescence to the that of the 12.5 ng/ml standard.



KIT COMPONENTS

Materials provided with the kit:

- Acridan & HRP conjugates, 5 **Store at -20°C**
- Dog PCT standard, 5 **Store at -20°C**
- Diluent: CSD50-1, 2 x 50 ml
- Background reducer: BR9-1, 9 ml
- Trigger solution: TS12-1, 12 ml
- 15 ml centrifuge tubes, 2

Materials required but not provided:

- VetBio-1 luminometer
- Precision pipettors and tips
- 12 x 75 mm borosilicate tubes²
- 12 x 75 mm test tube racks

STORAGE

Store the conjugate and PCT standard vials at or below -20°C. The remainder of the kit should be stored at 2-8°C.

GENERAL INSTRUCTIONS

1. Please take the time to completely read all instructions before starting your assay. Contact us if you need clarification.
2. All reagents used in the assay should be allowed to reach room temperature before use.

LUMINOMETER SETUP

1. Turn the VetBio-1 luminometer on.
2. Place the tubing from injector 1 into the tube containing background reducer. Tube holders are positioned adjacent to the injectors.
3. Place the tube from injector 2 into the tube containing trigger solution.
4. From the Protocol manager on the keypad select "Prime & Wash", "Prime", then "Start".
5. Open the luminometer drawer and insert an empty 12 x 75 mm tube.
6. Close the drawer and click "Start" again.
7. Wait for priming to complete, open the drawer and discard the tube.
8. From the Protocol manager select "Measure" followed by "Dog PCT".
9. Select "Start".
10. The experiment setup screen will be displayed, allowing you to enter experiment name, comments, and the number of samples. After entering the information select "Start".
11. Enter the sample IDs and dilution factor(s), if different than the default dilution factor.

¹ SPARCL technology, using acridan- and HRP-conjugated antibodies, was developed by and is licensed from Lumigen Corp.

² Only cylindrical 12 x 75 mm borosilicate glass tubes can be used. Do not use polystyrene or polycarbonate tubes.

12. The luminometer is now ready for use.
13. Press "Start" when you are ready to measure luminescence.

REPLICATES

The VetBio-1 luminometer protocols are configured to run singlets of the blank, standard and samples. At the discretion of the user it can be configured to run replicates. For the dog procalcitonin assay we recommend using duplicates.

1. From the protocol manager select "Measure" and the program you would like to modify.
2. Select "New" and "Copy Protocol".
3. Increase the Replicates as desired.
4. Create a new protocol name.
5. Select "Protocols" and save the new protocol.

SAMPLE PREP

This assay was designed for measurement of PCT in serum. In studies at Veterinary Biomarkers, Inc., we found PCT levels ranging from 0 to >100 ng/ml. Testing at a dilution of 20-fold is recommended. A 20-fold dilution can be obtained by mixing 12.5 µl of serum with 237.5 µl of diluent

STANDARD PREP

The PCT standard is provided ready to use. Thaw the vial shortly before use.

CONJUGATE PREP

1. The acridan and HRP conjugate mix should be prepared just before use (step 5 in the Procedure section).
2. Tap the vial to ensure that the contents are at the bottom of the vial before carefully removing the stopper.
3. Add 5.5 ml of diluent CSD50-1 to the vial. Insert the stopper and mix gently by inversion several times.
4. Each vial of reconstituted conjugate mix provides enough reagent to measure 0 and 12.5 ng/ml standards and up to eight samples (ten tubes) tested as singlets.³

PROCEDURE

1. Determine the number of 12 x 75 mm borosilicate glass tubes required for the assay. Ensure that all tubes fit easily in the VetBio-1 sample holder.
2. Pipet 100 µl of diluent into assay tube 1 (1 and 2, if testing duplicates). This serves as the blank.
3. Pipet 100 µl of the 12.5 ng/ml PCT standard into tube 2 (3 and 4, if testing duplicates).
4. Pipet 100 µl aliquots of the diluted samples into tubes 3, 4, 5... as defined by your assay format.
5. Add 0.50 ml of freshly prepared conjugate mix to each tube and mix gently. A vortex mixer may be used if available.
6. Incubate the mixtures at room temperature.
7. After 45 minutes insert tube 1 into the sample holder of the luminometer and close the drawer. The luminometer automatically injects background reducer and trigger solution, then measures luminescence (RLU/s).
8. Once the RLU/s value is recorded on the screen open the drawer and discard the tube.
9. Determine luminescence for the remaining tubes.
10. PCT concentrations are automatically calculated.
11. After measurement of the last sample, select "End".
12. Results will be saved but may be exported as Excel or pdf files via a USB stick.

LUMINOMETER MAINTENANCE

The luminometer injectors must be cleaned with distilled or deionized water at the end of each day the instrument is used to avoid clogging of the injector ports.

1. From the Protocol manager screen select "Prime & Wash".
2. Select "Backprime" to return unused reagents to the respective tubes.
3. If future use of the reagents is intended store the sealed tubes in a refrigerator.
4. Place the tubing from injectors 1 and 2 into separate 15 ml centrifuge tubes containing distilled or deionized water.
5. Select "Wash" from the "Prime & Wash" screen. Press "Start", insert an empty tube into the luminometer, close the drawer and press "Start".
6. Discard the tube.
7. Leave the injector tubing immersed in water.
8. Switch the luminometer off. It should not be left on.

³ If using multiple vials of conjugates to measure more than eight samples, combine the reconstituted contents of all vials and mix briefly before dispensing 0.5 ml aliquots into the reaction tubes.

ASSAY PERFORMANCE

Typical data: The table below shows results from an assay in which the zero standard (blank), 12.5 ng/ml standard, and three samples were tested in duplicate.



**Veterinary
Biomarkers, Inc.**

Dog PCT Exp#001


10.01.2020 22:04:39

Single Assay Start by : DOOR Blank : ON
 Measurement Delay [s] : 2.600 Time [s] : 0.800
 Injector 1 Delay [s] : 1.000 Volume [µL] : 50
 Injector 2 Delay [s] : 2.000 Volume [µL] : 100

Experiment comment: Test of three serum samples in duplicate

| Sample | Dilution | Rep | RLU/s | Conc [ng/mL] |
|----------|----------|-----|--------|--------------|
| Blank | | 1 | 2,612 | 0 |
| | | 2 | 3,230 | 0 |
| | | Avg | 2,921 | 0 |
| Standard | | 1 | 37,775 | 12.282 |
| | | 2 | 39,012 | 12.718 |
| | | Avg | 38,394 | 12.5 |
| 154 | 20 | 1 | 15,238 | 86.801 |
| | | 2 | 13,814 | 76.767 |
| | | Avg | 14,526 | 81.784 |
| 168 | 20 | 1 | 7,382 | 31.442 |
| | | 2 | 7,751 | 34.04 |
| | | Avg | 7,567 | 32.741 |
| 170 | 20 | 1 | 6,004 | 21.725 |
| | | 2 | 6,021 | 21.848 |
| | | Avg | 6,012 | 21.786 |

Linearity: To assess the linearity of the assay, serum samples containing PCT at concentrations of 38 and 100 ng/ml were serially diluted with diluent CSD50-1 to produce values within the dynamic range of the assay.



Rev 012020